FERNO Venice Carry Chair

4-wheeled carry chair



- **Lightweight and compact design** – Weighing just 9.9 kg, the Venice carry chair is easy to carry, with effortless manoeuvrability when using the chair. A simple to operate wall bracket is also available for storage.
- **Large rubberised wheels** – large high-grip rubberised wheels makes it easier to manoeuvre over a wide range of surfaces and environments, both in and outdoor
- **Folding footrest** – supporting the patient's feet and legs, with ankle strap for additional safety
- **Rear brakes for additional safety** – ensuring the chair is static and stable when transferring someone on or off the chair

The enhanced model of carry chair from Ferno, the Venice chair.

The Venice carry chair is much more functional than other carry chairs and is designed to support a greater physical size and weight range, with a SWL of 230 kg (36 stone). Weighing in at only 9.9 kg, the chair is also easy to carry. Combined with the 4 rubberised wheels, 2 of which rotate at the front, it is easy to operate and has excellent manoeuvrability and stability. Brakes on the rear wheel ensures full safety when transferring someone on or off the chair. Folding lift handles at the rear of the chair aid when wheeling the chair as well as making it easier to lift up over curbs, thresholds and other obstacles, with less effort. The chair neatly and easily folds for storage, measuring only 20 cm wide once folded up.





FERNO Venice Carry Chair







Features & Benefits

- **4-wheeled chair** – large rubberised profiled wheels with enhanced manoeuvrability
- **Rear fold-out manoeuvre handles** – angled for improved ergonomics when wheeling the chair or needing to lift over a threshold or other obstacle
- **Ergonomic Frame -** lightweight chair frame easy to lift or carry and with SWL of 230 kg (36.2 stone)
- **Improved seat design -** seat and backrest slightly angled back for better positioning and comfort for the patient
- **Lower-height footrest** – easier to place feet, especially useful for bariatric patients
- **Easy and compact to store** – no need to lift to place in bracket for storage, folds to a depth of 20 cm, and with magnets to hold the wheels in place

Optional Accessories

- **Options** – various optional extras also available according needs of use, including backrest and head cushion, extended leg support

Specifications

HEIGHT - Deployed	859 mm
HEIGHT -Stowed	908 mm
WIDTH	503 mm
DEPTH - Handles Deployed	990 mm
DEPTH - Handles folded	763 mm
DEPTH - Chair folded	203 mm
SEAT HEIGHT	484 mm
SWL	230 kg / 36.2 stone
CHAIR WEIGHT	10 kg
SimpleClack storage bracket	3.4 kg



- Mesh back **support** giving support to the head or shoulders for taller patients



- Padded head **support** giving additional sup-port/comfort for patient's head



- 'Simple Clack' Storage Bracket securing the chair safely against a wall



Extended Leg Support Option - transport patient with injured leg

Part Numbers

21-0081-001	Venice carry chair, with webbing straps + backrest cushion
25-1000-020-UK	SimpleClack in-vehicle/wall storage bracket
10-1964-001	Mesh netting back/head support
031385825	Padded head support
25-00285	Extended leg support option (order at time of purchasing chair)

Supplied Backrest cushion



FERNO

Call +44 (0) 1274 851 999 **email** sales.uk@ferno.com **www** ferno.com

Ferno (UK) Limited, Stubs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ, England